



## Visionary Natural Starsoft 66 UV / Starsoft 75 UV

### Geometrie:

Spherical front and back optic Zones  
Aspheric flattening  
Optimized Opticzone



### Technical data:

Total diameter: 13,00 -> 18,00 mm  
Base curve: 8,00 -> 10,00 mm  
Sphere: -40,00 -> +40,00 dpt  
Steps 0,10 mm

### Contact Lens Care Instructions

Important: Especially contact lens care products that contain EDTA will quickly destroy the coloration of your Visionary lens. Please use one step-T solution

Irisdiameter: 10,5/ 11,0/ 11,5/ 12,0/ 12,5  
Clear pupil: 3,0/ 3,5/ 4,0/ 4,5/ 5,0/ 5,5/ 6,0  
Black pupil: 3,0/3,5/ 4,0/ 4,5/ 5,0/ 6,0



Star-Lens AB • Box 170 • SE-662 24 Åmål  
Tel +46 (0)532 607980 • Fax +46 (0)532 185 25 • [bestallning@starlens.se](mailto:bestallning@starlens.se)  
[www.starlens.se](http://www.starlens.se)